

pycaspule - Python package

Impulse Neiry

May 5, 2021

Contents

1	Introduction	3
1.1	Installation	3
2	Reading data with RecordReader	3
3	Importing data into MNE	5

1 Introduction

We provide a `pycapsule` package with various tools for working with Capsule BCI system. At this moment it allows to:

1. Read session records into memory
2. Import data into MNE

1.1 Installation

Python package can be installed via pip:

```
pip install pycapsule
```

Alternatively it can be found at `{SDK-ROOT}/Tools/Python/Packages/pycapsule`.

2 Reading data with RecordReader

Record reader functionality can be used to unpack stream data from sessions ran by Capsule. Reader works with session data files (`.dat/.rec`) in pairs. Both are required in order to successfully unpack data. Unpacking data using our package requires implementing a subclass of `RecordReaderVisitor`. It is a class which will receive all incoming packets as they are unpacked.

- `P300ProcessingUnit` packet contains all the data necessary to perform model training and predictions.
- `Raw EEG` and `Raw Resistance` packets contain continuous signal data for the whole session.
- `Interface` data packets contain specific commands related to a particular interface, these commands are issued by client.

Let's look at a simple example of how to unpack session data:

Listing 1: Usage example

```
1 from pycapsule.record_reader import *
2
3 class Visitor(RecordReaderVisitor):
4     def OnP300ProcessingUnit(self, p300unit:P300ProcessingUnit):
5         print("P300 Unit!")
6     def OnRawEEG(self, eegData:np.ndarray, eegTimestamps:np.array):
7         print("Raw EEG!")
8     def OnRawResistance(self, resData:np.ndarray):
9         print("Raw Resistance!")
10    def OnInterfaceData(self, interfaceData):
11        print("Interface Data!")
12
13 visitor = Visitor()
14
15 filepath = "D:/Capsule-recs/Impulse - 2020.12.07 - 07.56.32/session.rec"
16 RecordReader.Unpack(filepath, visitor)
```

We will now review steps required to retrieve data:

- line 1 Import Capsule's record_reader required for data unpacking
- line 3 Implement subclass of a RecordReaderVisitor, there is no need need to implement every function
- line 16 Call unpack functionality with a filepath to a session recording (.rec data file) and instance of a visitor we implemented

Here's a basic reference of what data classes contain:

Listing 2: Data classes

```
1 @dataclass
2 class P300StimulusData:
3     stimulusId: int
4     sampleStartIndex: int          # index into eegData
5     timestamp: float              # in seconds
6
7 @dataclass
8 class P300ProcessingUnit:
9     unitId: int
10    actId: int
11    targetStimulus: int
12    stimuliCount: int
13    eegData: np.ndarray             # shape(channels, samples)
14    eegTimestamps: np.array         # shape(samples)
15    stimuliData: typing.List[P300StimulusData]
16    shouldEndLearn: bool           # (legacy) signals BCI system to train a model
```

P300ProcessingUnit in a basic sense represents a sequence of stimuli activations. Stimuli with their id, timestamp and index into eegData are stored in stimuliData. EEG data for a sequence is stored in eegData numpy array with shape (channels, samples).

3 Importing data into MNE

It is also possible to import records directly into MNE.

Listing 3: Usage example

```
1 from pycapsule.mne.io import read_raw_csr
2 import matplotlib
3
4 filepath = "D:/Capsule-recs/Impulse - 2020.12.07 - 07.56.32/session.rec"
5 raw, events, event_id = read_raw_csr(filepath)
6 raw.set_montage('standard_1020', match_case=False) # montage to visualize electrode placement
7 raw.filter(1, 40) # since data is not filtered, we apply basic filter
8 raw.plot(events, event_id=event_id, event_color={1: 'g', 2: 'b', 3: 'r'})
9 raw.plot_psd()
10 matplotlib.pyplot.show()
```




Figure 1: Plotted filtered signal on O1 channel with markers




Figure 2: PSD plot with montage